

DVMAXX HD

Product Features

The DVMAXX HD is the latest high-definition medical video recorder and image capture designed for quick set up and ease of use. With auto sensing video connections and automatic video resolution detection, set up is quick and easy. The front control panel allows users to intuitively record videos and photos in crystal clear HD, and a built in hard drive allows for simultaneous recording to protect your data.



- Full HD 1920 x 1080 video recording
- 3.5" LCD monitor for live view and playback
- Built in hard drive for simultaneous recording to USB
- Network capable
- Automatically detects connections and video resolution

SPECIFICATIONS

Recordable Devices	USB Flash Drive, USB Hard Drive, Internal Hard Drive (500GB)
Display	3.5 inch LCD
Power Supply	AC 100V - 240V 50/60Hz
Operating Temperature	41 to 104 deg F (5 to 40 deg C)
Operating Humidity	5 to 80% (No condensation)
Dimension	210.5mm W x 235mm D x 88.5mm H, 8.25in W x 9.25in D x 3.375in H
Weight	2.6kg, 5.7lbs
Image Format	Jpeg
Video Format	MPEG-4 AVC/H.264
Video Input	SDI, DVI, S-Video, Composite (BNC)
Video Output	SDI, DVI, S-Video, Composite (BNC), SD/HD-SDI
Audio Input	SDI, Stereo mini jack
Audio Output	SDI, Stereo mini jack
Other Interfaces	RJ-45, USB 2.0 Type A, UBS 2.0 Type B, Stereo mini jack
Regulatory Compliance	IEC60601-1, IEC60601-1-2, CSA60601-1, EN60601-1, EN60601-1-2, ICES-003, VC CI Class B, FCC Class B, FDA Class 1, ISO13485, CSA NRTL, CE