


MEDVIX™ - AMVX2608HD

High quality surgical display with HD Picture.

Medvix AMVX2608HD is a 26-inch LED backlight Multi-Modality display that is purposely built for the operating room. The AMVX2608HD offers an easy to clean design, smart mechanics and it provides the best detailed images for any surgical suite.

- 26" IPS-Pro wide-screen LCD
- High-bright LED backlight with output stabilization
- Full HD resolution (1920 × 1080)
- Wide viewing angle thanks to IPS-Pro LCD technology
- 10-bit image processing
- 12-bit look-up table
- Artifact-free images
- Calibrated, consistent colors and grayscales
- Easy-clean design with sealed front (IPX5)
- Smart cable management system



PICTURE PERFORMANCE

LCD Panel	26" TFT AM LCD / IPS-Pro technology / LED backlight
Resolution	2MP (1920 x 1080)
Pixel Pitch	0.3 mm
Luminance	450 cd/m ²
Aspect Ratio	16:9
Contrast Ratio	1400:1
Viewing Angle (H/V)	178° H / 178° V

INPUT/OUTPUT

Video Inputs	2X DVI-I Single link (with HDMI support), RGBS / YPbPr, S-video, Composite Video, 3G-SDI, VGA (on DVI-I),
Video Outputs	DVI-D, S-video, Composite Video, 3G-SDI

GENERAL

Remote Control	RS-232
Power Requirements	100–250 V / 47–63 Hz auto-switch
Mounting Standard	VESA (100 x 100 mm)
Dimensions (W × H × D)	633 × 411 × 90 mm (24.9 × 16.2 × 3.5")
Safety / EMC	IEC60601-1, UL60601-1, CSA-C22.2, Protection level: IP21 (with Protection screen) / Front face: IPX5 protection level Green requirements: ROHS2, REACH, WEEE compliant
Weight	9.5 kg (20.9 lbs) (Monitor only)
OSD Languages	English, French, German, Spanish, Italian
Recommended Modalities	Endoscopy, US, X-ray, Surgery, CR, MR, CT, NM
Warranty	2 years



©2018 AMPRONIX INC. 14002