



MGC

LAMPS

Part of MGC Lighting Group

Medical Batteries





YEARS OF LIGHTING

Philips

Yuasa

Maquet

Welch Allyn

**Duracell
Industrial**

Panasonic

Multipower



MGC

LAMPS

MGC Lamps is proud to celebrate 40 years of professional and unrivalled service.

We have come a long way in those 40 years, from humble beginnings supplying film & television studios, to being one of the largest independent lighting distributors in Europe.

The last 20 years has seen MGC Lamps (Lighting Group) emerge as one of the largest lamp suppliers to the NHS, supplying everything from xenon lamps, lamps for dermatology, lamps for operating theatres, hand held equipment, and much more.

We don't only do lamps!

MGC has expanded its portfolio over the years and now stocks associated equipment and products, from Keeler to fibre optic cables. MGC can offer a wide range of products at highly competitive pricing.

Our range of batteries is testament to the growth and commitment of MGC in medical and scientific products.











If you cannot find what you need, please ask.

Index

Manufacturer	Page
AED	2
Aesculap	2
Aestiva	2
Alaris	2
Arjo	2-3
Baxter	3
Braun	3
Cardiac Science	3
Cardiette	3
Cardioline	3-4
Covidien	4
Criticare	4
CU Medical Systems	4
Datascope	4-5
Datex Ohmeda	5-6
Draeger	6
Fiamm	6
Fresenius	7
Fukuda Denshi	7
GE Criticon	7
GE Marquette	7-9
Graseby	9
Heine	9
Hewlett Packard	9-10
Ivac	10-11
Laerdal	11
Maquet	11
Masimo	11
Nellcor	11-12
Nihon Kohden	12
Nonin	12
Philips	12-13
Physio Control	13-14
Protocol	14
Welch Allyn	14-15
Zoll	15
Duracell Industrial Batteries	16
Terms and Conditions	17













The indicated brands and company designations are in parts, registered trademarks of their respective owners and are only used to indicate their compatibility.

Manufacturer	System	Model	Spec	Property	Original or Non Original	MGC Part No.	
AED	AED Defib	Defibtech Lifeline View Pro, ReviveR, ECG Type DCF-2003	12V 2.8Ah	Lithium	Original	BATT/110624-O	no image available
Aesculap	Electro Stimulator		6V 0.18Ah	NC	Non original	BATT/110181	
Aestiva	Aestiva Anaesthetic Machine	Yuasa NP2.3-12	12V 2.3Ah	Pb	Non original	BATT/361023	
Alaris	Alaris Medical System (Carefusion) Asena Syringe Pump	Type CC/GH/PK/GS, Smart Pack 1000SP01122 1000SP01080 1000SP01798 (with PCB)	7.2V 2.7Ah	NiMH	Non original	BATT/110143-K	
Alaris	Alaris Medical System (Carefusion) Asena Infusion Pump	Type 1000SP00487, 1000SP01302, 1000SP01794	9.6V 2.7Ah	NiMH	Non original	BATT/110324	
Alaris	Alaris Medical System Infusion Pump Asena GW		7.2V 1.25Ah	NiMH	Non original	BATT/110178	
Alaris	Alaris Infusion Pump TIVA		6V 3.0Ah	Pb	Non original	BATT/Z-366030-1	
Alaris	Alaris Medical System Gateway Workstation	Type 1000SP00605	9.6V 2.7Ah	NiMH	Non original Retrofit	BATT/110475-U	
Arjo	Arjo Lifter	Alenti/Encom/Miranti/Marisa Type NDA0100-03	24V 4.0Ah	Pb	Non original Retrofit	BATT/110179-A	
Arjo	Arjo Lifter	Sara 3000 MaX/Move Type NDA0100-20	24V 4.0Ah	Pb	Non original Retrofit	BATT/110212-A	
Arjo	Arjo Lifter	Alenti/Encom/Miranti/Marisa Sara 3000/Bolero Type NDA0100-03, NDA0100-20, NDA0100-32	24V 4.0Ah	Pb	Original	BATT/110212-O	

01473 466 300



MGC Lamps Ltd
 1 Sovereign Centre,
 Farthing Road, Ipswich
 Suffolk, IP1 5AP

Manufacturer	System	Model	Spec	Property	Original or Non Original	MGC Part No.	
Arjo	Arjo Lifter	Trixie/Sarita Type KTA0102	24V 2.0Ah	Pb	Non original Retrofit	BATT/110237-A	
Arjo	Arjo Lifter	Trixie/Sarita Type KTA0102	24V 2.0Ah	Pb	Original	BATT/110237-O	
Arjo	Arjo Lifter	Tempo/Opera Type KPA0100	24V 5.0Ah	Pb	Non original Retrofit	BATT/110325-A	
Baxter	Baxter Pulsoximeter	ASAT	8V 2.5Ah	Pb	Non original	BATT/110283	
Braun	Braun Perfusor	FT/V11	7.2V 1.5Ah	NC	Non original	BATT/110009	
Braun	Braun Perfusor	FT91/F/Secura P	7.2V 1.2Ah	NC	Non original	BATT/110010	
Braun	Braun Perfusor	FM (MFC)	7.2V 1.2Ah	NC	Non original	BATT/110013	
Braun	Braun Perfusor	Infusomat Space, 8713180, 34520856	4.8V 2.1Ah	NiMH	Original	BATT/110256	
Cardiac Science	Cardiac Science Powerheart AED G3 Defib	Type 9145/9146/9300E/9300P	12V 7.5Ah	Lithium	Original	BATT/110313	
Cardiette	Cardiette ECG Recorder	AR600ADV	9.6V 2.4Ah	NiMH	Non original	BATT/110236	
Cardioline	Cardioline Cardiette ECG Recorder	AR1200, AR1200ADV, AR1200view, FC1400	12V 1.5Ah	NiMH	Non original	BATT/110176	

United Kingdom Sales

☎ 01473 466 300
☎ FREEPHONE 0800 137 739
✉ enquiry@mgc-lighting.com

International Sales

☎ +44(0) 1473 462321
✉ internationalsales@mgc-lighting.com




Can't find what you need?
Just call!

Manufacturer	System	Model	Spec	Property	Original or Non Original	MGC Part No.	
Cardioline	Cardioline ECG	Delta, Delta1+, Delta3+	24V 0.7Ah	NC	Non original	BATT/110222	
Cardioline	Cardioline ECG	Delta 60, 60 Plus	24V 1.3Ah	NC	Non original	BATT/110260	
Covidien	Covidien Aerogen Aereoneb Pro	Type AP0200	4.8V 0.45Ah	NiMH	Non original	BATT/110349	
Criticare	Criticare Pulsoximeter	502, 504, 504US, 50245, 50445, 506, DXN	6V 4.5Ah	Pb	Non original	BATT/301010	
CU Medical Systems	Paramedic Defib	CU-ER Series type CU-ER1	12V 4.5Ah	NiMH	Original	BATT/110615-01	
Datascope	Datascope Mindray Monitor T5/T8	Type M05-010002-06	11.1V 4.5Ah	LiOn	Non original	BATT/110417	
Datascope	Datascope Mindray Monitor T5/T8	Type M05-010002-06	11.1V 4.5Ah	LiOn	Original	BATT/110417-0	
Datascope	Datascope Accutorr Plus/Trio Patient Monitor	0146-00-0069 (0146-00-0099) 115-018015-00 LI23S003A	11.1v 4.4Ah	LiOn	Original	BATT/110262	
Datascope	Datascope Monitor Duo	0146-00-0079	7.4V 7.8Ah	LiOn	Original	BATT/110339	
Datascope	Datascope Mindray Vital Functions Monitor	VS800, PM-8000, PM-9000, MEC-1200VET, IPM9800	11.1v 4.4Ah	LiOn	Original	BATT/110404-0	
Datascope	Datascope Mindray PM60 Pulsoximeter L/115001A	Type M05-010004-08	3.7V 1.8Ah	LiOn	Original	BATT/110492-0	

01473 466 300



MGC Lamps Ltd
 1 Sovereign Centre,
 Farthing Road, Ipswich
 Suffolk, IP1 5AP

Manufacturer	System	Model	Spec	Property	Original or Non Original	MGC Part No.	
Datascope	Datascope Mindray Monitor	PM8000 Express/PM9000/MEC1000/MEC1299 Type 49471-4	12V 2.3Ah	Pb	Non original	BATT/399013	
Datex Ohmeda	Datex Ohmeda/Braun Pulsoximeter	Biox 3700E	8V 2.5Ah	Pb	Non original	BATT/110058	
Datex Ohmeda	Datex Ohmeda/Braun Pulsoximeter	3800	8V 3.0Ah	Pb	Non original	BATT/110059	
Datex Ohmeda	GE/Datex Ohmeda Monitor	F/FM/F-FM-01/ F-CUB.07 Series/ Type M10008142	11.1V 3.6Ah	LiOn	Non original	BATT/110481	no image available
Datex Ohmeda	GE/Datex Ohmeda Monitor	F/FM/F-FM-01/ F-CUB.07 Series/ Type M10008142	11.1V 3.52Ah	LiOn	Original	BATT/110481-O	
Datex Ohmeda	Datex Ohmeda Pulsoximeter	Biox 3770/3775	9.6V 1.4Ah	NC	Non original	BATT/110109	
Datex Ohmeda	Datex Monitor (takes 2 batteries)	S5, Type 17014	12V 4.0Ah	NiMH	Non original	BATT/110141	
Datex Ohmeda	Datex Ohmeda Monitor	Cardiocap 5	12V 2.6Ah	Pb	Non original	BATT/110240	
Datex Ohmeda	Datex Ohmeda Light Monitor	Type 893365	13.2V 1.1Ah	NC	Non original	BATT/110269	
Datex Ohmeda	Datex Ohmeda Light Monitor	Type 896895, 887753	13.2V 7.2Ah	NC	Non original Insert only - no case	BATT/110197	
Datex Ohmeda	Datex Ohmeda TruSat Pulsoximeter	6050-0006-578 REV J	4.8V 3.2Ah	NiMH	Non original	BATT/110335	











United Kingdom Sales

☎ 01473 466 300
☎ FREEPHONE 0800 137 739
✉ enquiry@mgc-lighting.com

International Sales

☎ +44(0) 1473 462321
✉ internationalsales@mgc-lighting.com

Can't find what you need?
Just call!

Manufacturer	System	Model	Spec	Property	Original or Non Original	MGC Part No.	
Datex Ohmeda	GE Datex Ohmeda Monitor	B30	11.1V 3.6Ah	LiOn	Non original	BATT/110487	
Datex Ohmeda	GE Datex Ohmeda Monitor	B30	11.1V 3.52Ah	LiOn	Non original	BATT/110487-0	
Draeger	Siemens/Draeger Monitor Backup Battery	Infinity Gamma/Gamma XL (Backup battery)	3V	Lithium	Non original	BATT/110357	
Draeger	Draeger Monitor	PM 8060, Type 1835343	3V	Lithium	Non original	BATT/110483	
Draeger	Draeger Monitor	PM 8060, Type 1835343	3V	Lithium	Original	BATT/110483-0	
Draeger	Draeger Ventilator	Oxylog 2000, 8411599	7.2V 1.0Ah	NC	Non original	BATT/110146	
Draeger	Draeger Ventilator	Oxylog 2000 PLUS, 3000, 3000 PLUS Type 2M86733	11.1V 7.2Ah	LiOn	Original	BATT/110242-0	
Draeger	Draeger Ventilator	Draeger Primus Evita 4 Panasonic LC-XD1217P Screw connection	12V 17.0Ah	Pb	Original Panasonic	BATT/301003	
Draeger	Draeger Ventilator	Oxylog 2000 PLUS, 3000, 3000 PLUS Type 2M86733	10.8V 8.7Ah	LiOn	Non original	BATT/110242	
Draeger	Draeger Monitor	Infinity Gamma/ Gamma XL Patient Monitor (new version 2014)	14.4V 5.9Ah	LiOn	Original	BATT/110196-01	no image available
Fiamm	Fiamm Filtration System	Filtration System used in CCU FG20121	12V 1.2Ah	Pb	Original	BATT/371012	

01473 466 300



MGC Lamps Ltd
 1 Sovereign Centre,
 Farthing Road, Ipswich
 Suffolk, IP1 5AP

Manufacturer	System	Model	Spec	Property	Original or Non Original	MGC Part No.	
Fresenius	Fresenius Infusion Pump	Vial/MCM/ Volumed VP2000/5000	12V 1.70Ah	NC	Non original	BATT/110022	
Fresenius	Fresenius Infusion Pump	Volumat Agilla type 179033-R2	7.2V 2.2Ah	LiOn	Non original	BATT/110320	
Fukuda Denshi	Fukuda Denshi ECG	ECG CardiMax FX-7202, Type T8HRAAU-4713	9.6V 1.75Ah	NiMH	Non original Retrofit	BATT/110490	
GE Criticon	GE Carescape VSM	Procure Dinamap 100/300/400/420/V100 GH633/633178	6V 3.3Ah	Pb	Non original	BATT/110564	
GE Marquette	GE Hellige ECG Recorder	MAC PC, MAC PC6, MAC 12	19.2V 1.70Ah	NC	Non original	BATT/110150	
GE Marquette	GE Hellige ECG Recorder	MAC 8, MAC Personal	19.2V 4.50Ah	NC	Non original	BATT/110162	
GE Marquette	GE Hellige Marquette Patient Monitor	MAC3500/MAC5000ST/MAC5500/Dash 3500/Dash 5500 Type 900770-001	18V 3.5Ah	NiMH	Non original	BATT/110193	
GE Marquette	GE Hellige Marquette ECG Cardiosmart	MAC 1200 ST, 1100, PC1000	18V 1.3Ah	NC	Non original	BATT/110032-I	
GE Marquette	GE Hellige Marquette Patient Monitor		13.2V 0.8Ah	NC	Non original	BATT/110274	
GE Marquette	GE Hellige Marquette Eagle Monitor	1000/1006/1008/1009	12V 2Ah	NC	Non original	BATT/110107	
GE Marquette	GE Hellige Marquette Eagle Monitor	4000	12V 2.5Ah	NC	Non original	BATT/110184	

United Kingdom Sales

☎ 01473 466 300
☎ FREEPHONE 0800 137 739
✉ enquiry@mgc-lighting.com

International Sales

☎ +44(0) 1473 462321
✉ internationalsales@mgc-lighting.com

Can't find what you need?
Just call!

Manufacturer	System	Model	Spec	Property	Original or Non Original	MGC Part No.	
GE Marquette	GE Marquette Defib	900/1200/1250/1500	24V 1.70Ah	NC	Non original Retrofit	BATT/110200	
GE Marquette	GE Hellige Marquette Monitor	Dash 2000	12V 3.0Ah	NC	Non original	BATT/110202	
GE Marquette	GE Hellige Marquette Monitor	Dash 3000/4000/5000 Transport Pro B30 SM201-6. Type 2017857-002	11.1V 3.6Ah	LiOn	Non original	BATT/110155-SM	
GE Marquette	GE Hellige Marquette Monitor	Dash 3000/4000/5000 Transport Pro B30	11.1V 3.6Ah	LiOn	Original	BATT/110155-O	
GE Marquette	GE Marquette Monitor	Carescape B650 2036984-001	11.1V 6.0Ah	LiOn	Original	BATT/110465-O	
GE Marquette	GE Marquette Monitor	Carescape Patient Data Module 2016989-002	10.8V 1.8Ah	LiOn	Original	BAT/110402-O	
GE Marquette	GE Marquette Monitor	Dash 2500 Type 2023852-029/ 2023227-001	8.4V 8.0Ah	NiMH	Original	BATT/110450-O	
GE Marquette	GE Marquette Hellige Defib	Responder 3000/Dash Responder/ Type 30344631/ Part 2009218-001	12V 2.0Ah	NC	Original	BATT/110131	
GE Marquette	GE Marquette Hellige Defib	Responder 2000 Type 2025267-001	11.1V 6.6Ah	LiOn	Original	BATT/110396	
GE Marquette	GE Marquette Defib	Responder AED 2019437-001		Lithium	Original	BATT/110329	
GE Marquette	GE Marquette Hellige Operating Table		12V 8.0Ah	Pb	Non original	BATT/110130	

01473 466 300



MGC Lamps Ltd
 1 Sovereign Centre,
 Farthing Road, Ipswich
 Suffolk, IP1 5AP

Manufacturer	System	Model	Spec	Property	Original or Non Original	MGC Part No.	
GE Marquette	GE Marquette Hellige Operating Table	Type 1150.01.A1	12V 7.0Ah	Pb	Non original	BATT/110130-2	
Graseby	Graseby Infusion Pump	Cyclon D-Cell, 2000/3000 Series, 0810-0004	2V 2.5Ah	Pb	Non original	BATT/399009	
Heine	Heine Ophthalmoscope	S3Z, X-002.99.314	3.6V 0.65Ah	NiMH	Original	BATT/110901	
Heine	Heine Ophthalmoscope	S3Z, X-002.99.314	3.6V 0.65Ah	NiMH	Non original	BATT/110901-AN	
Heine	Heine Ophthalmoscope	Beta X-02.99.382, X-02.99.380	3.5V 1.0Ah	NiMH	Original	BATT/110904	
Hewlett Packard	Hewlett Packard Defib	Codemaster 100 Type M2476B	12V 3.0Ah	NC	Non original Retrofit	BATT/110097-I	
Hewlett Packard	Hewlett Packard Defib	Codemaster 100 Type M2177B/M2477B	12V 4.0Ah	NC	Non original Retrofit	BATT/110098-I	
Hewlett Packard	Hewlett Packard Defib Monitor	Codemaster XL Plus M1722A/B, MN1723A/A, 1724A M1758A	12V 4.2Ah	Pb	Non original	BATT/110039	
Hewlett Packard	HP/Philips Pagewriter ECG Machine	100, 200, 300 i. 300Pi, M1770A/M1771A/M1772A	6V 7.2Ah	Pb	Non original	BATT/110244	
Hewlett Packard	HP Transport Monitor	40488A/M1275, Flat battery - connection on the side	12V 2.10Ah	Pb	Non original	BATT/399013	
Hewlett Packard	Hewlett Packard Monitor	M3046 Vindia M3/M4 Type NJ1020HP	12V 3.8Ah	NiMH	Non original	BATT/110001-X	











United Kingdom Sales

☎ 01473 466 300
☎ FREEPHONE 0800 137 739
✉ enquiry@mgc-lighting.com

International Sales

☎ +44(0) 1473 462321
✉ internationalsales@mgc-lighting.com

Can't find what you need?
Just call!

Manufacturer	System	Model	Spec	Property	Original or Non Original	MGC Part No.	
Hewlett Packard	Hewlett Packard Defib/ Monitoe	43100A/43110A/43120A/43130Ah/43200A 123	12V 2.5Ah	Pb	Non original	BATT/110038	
Hewlett Packard	Hewlett Packard Monitor	78660A-60401, 78670, 671, 672	12V 2.3Ah	NC	Non original	BATT/110623	no image available
Ivac	Ivac Alaris Syringe Pump	741, Panasonic LC-R061R3P	6V 1.3Ah	Pb	Original Panasonic	BATT/301006	
Ivac	Ivac Alaris Syringe Pump	P1000-P7000	6V 3.0Ah	Pb	Non original	BATT/366030	
Ivac	Ivac Alaris Infusion Pump	Signature Gold 7101/7130/7131/7230/7231	12V 1.8Ah	NC	Non original	BATT/110113	
Ivac	Ivac Infusion Pump	231/531/281	12V 0.6Ah	NC	Non original	BATT/110114	
Ivac	Ivac Infusion Pump	561/565/570/571/572	12V 2.5Ah	Pb	Non original	BATT/110043	
Ivac	Ivac Infusion Pump	581/591/597/598	6V 3.4Ah	Pb	Non original	BATT/110044	
Ivac	Ivac Syringe Pump	707, 711	8V 2.5Ah	Pb	Non original	BATT/110045	
Ivac	Ivac Syringe Pump	770	6V 2.5Ah	Pb	Non original	BATT/110046	
Ivac	Ivac Syringe Pump	770 (Backup)	4.8V 0.16Ah	NC	Non original	BATT/110047	

01473 466 300



MGC Lamps Ltd
 1 Sovereign Centre,
 Farthing Road, Ipswich
 Suffolk, IP1 5AP

Manufacturer	System	Model	Spec	Property	Original or Non Original	MGC Part No.	
Ivac	Ivac Thermometer	Thermometer 2000	4.8V 0.7Ah	NC	Non original	BATT/110048	
Laerdal	Laerdal Suction Pump	790009, 790010, 790012, 790013, 790017 (takes 2 batteries)	6V 1.40Ah	NC	Non original	BATT/110082	
Laerdal	Laerdal Suction Pump	LCSU4 Type 886113	12V 1.76Ah	NiMH	Original	BATT/110517	
Laerdal	Laerdal Suction Pump	LSU/LSU4000 Type 780000, 780010, 780020	12V 2.0Ah	Pb	Non original	BATT/399035	
Laerdal	Laerdal Suction Pump	LSU 780800	12V 2.0Ah	NiMH	Original	BATT/780800	no image available
Maquet	Maquet Servo Ventilator System	Servo I, Servo S, 6487180	12V 4.0Ah	NiMH	Original	BATT/110203-0	
Maquet	Maquet Servo Ventilator System	Servo I, Servo S, 6487180	12V 3.8Ah	NiMH	Non original	BATT/110203	
Masimo	Masimo Pulse Oximeter	Radical 7 Hardpack	4.8V 1.8Ah	NiMH	Original	BATT/110345-0	
Masimo	Masimo Pulse Oximeter	Radical 7 Softpack	4.8V 1.6Ah	NiMH	Original	BATT/110409-0	
Nellcor	Nellcor Patient Monitor	N5500, N5600	14.4V 4.0Ah	NiMH	Non original	BATT/110273	

United Kingdom Sales

☎ 01473 466 300
 ☎ FREEPHONE 0800 137 739
 ✉ enquiry@mgc-lighting.com

International Sales

☎ +44(0) 1473 462321
 ✉ internationalsales@mgc-lighting.com

Can't find what you need?
Just call!

Manufacturer	System	Model	Spec	Property	Original or Non Original	MGC Part No.	
Nellcor	Nellcor Pulsoximeter	N10	6V 0.6Ah	NC	Non original	BATT/110194	
Nellcor	Nellcor Pulsoximeter	N560 Oximax, Mediana M6008-O	9.6V 3.8Ah	NiMH	Non original	BATT/110309	
Nihon Kohden	Nihon Kohden Patient Monitor	Lifescop 1. BSM2300, 2301A, 2301C, 2301-K, 2303C, 2303K, 2304A, 2351A, 2351C, 2351K, 2353C, 2353K, 2354A type X062, YS-076P5	12V 4.5Ah	NiMH	Non original	BATT/110216	
Nonin	Nonin Advant Pulsoximeter	2120/9600 type 4032-001	7.2V 4.0Ah	NiMH	Non original	BATT/110268	
Philips	Philips/Laerdal Monitor/Defib	HeartStart MRX	14.4V 6.3Ah	LiOn	Original	BATT/110245	
Philips	Philips/Laerdal Defib	HeartStart XL/XLT/4000/ M3535A/ M3536A/ M3500B/ M5500B Type M3516A/ M5516A	12V 2Ah	Pb	Non original	BATT/110271-A	
Philips	Philips/Laerdal Defib	HeartStart XL/XLT/4000/ M3535A/ M3536A/ M3500B/ M5500B Type M3516A/ M5516A	12V 2Ah	Pb	Original	BATT/110271-O	
Philips	HP/Philips Cardiograph	Pagewriter 100/200/300/300PI/ M1770A/M1771A/M1772A Type M2460A	6V 7.2Ah	Pb	Non original	BATT/110244	
Philips	HP/Philips Cardiograph	Pagewriter Trim I/Trim II/ Trim III Type PH98980313015	12V 2.6Ah	Pb	Non original Retrofit	BATT/110297-U	
Philips	HP/Philips Cardiograph	Pagewriter Trim I/Trim II/ Trim III Type PH98980313015	12V 2.9Ah	Pb	Original	BATT/110297-O	
Philips	HP/Philips Cardiograph	Pagewriter 10/10i type M2662A/M3941A	7.2V 2.7Ah	NiMH	Non original	BATT/110395	

01473 466 300



MGC Lamps Ltd
 1 Sovereign Centre,
 Farthing Road, Ipswich
 Suffolk, IP1 5AP

Manufacturer	System	Model	Spec	Property	Original or Non Original	MGC Part No.	
Philips	Philips Monitor	A3/3928A Type M3926	6V 8.0Ah	Pb	Non original	BATT/110275	
Philips	Philips Monitor	C3 Type86247 4C3/832478/453563480031	12V 4.50Ah	Pb	Non original	BATT/110276	
Philips	Philips Monitor	MP20/MP30/MP40/MP5/ MP50/MP70 M8001A/ M8002A. Type M4605A	10.8V 8.7Ah	LiOn	Non original	BATT/110219	
Philips	Philips Monitor	MP20/MP30/MP40/MP5/ MP50/MP70 M8001A/ M8002A. Type M4605A	10.8V 6.0Ah	LiOn	Original	BATT/110219-0	
Philips	Philips/laerdal Defib	Heartstart Forerunner 2/2* /FR2/FR2* M3840A/ M3841A/M3860A/M3861A Type M3863A	12V 4.2Ah	Lithium	Non original	BATT/110218	
Philips	Philips/laerdal Defib	Heartstart Forerunner 2/2* /FR2/FR2* M3840A/ M3841A/M3860A/M3861A Type M3863A	12V 4.2Ah	Lithium	Original	BATT/110218-0	
Philips	Philips Defib	Heartstart HS1, FRx non-rechargeable M5070A	9V 4.2Ah	Lithium	Original	BATT/110300-0	
Philips	Philips SureSigns Monitor	SureSigns VM3, 4,6,8, VS3. Type 989803194541, 989803144631	11.1V 7.2Ah	LiOn	Original	BATT/110348-0	
Physio Control	Physio Control Defib	Lifepak 20	12V 3.0Ah	NiMH	Original	BATT/110220	
Physio Control	Physio Control Defib	Lifepak 20	12V 3.0Ah	NiMH	Non original	BATT/110220-A	
Physio Control	Physio Control Defib	Lifepak 20e 3205296-002 Type 11141-000112	11.1V 6.0Ah	LiOn	Original	BATT/110220-L	











United Kingdom Sales

☎ 01473 466 300
☎ FREEPHONE 0800 137 739
✉ enquiry@mgc-lighting.com

International Sales

☎ +44(0) 1473 462321
✉ internationalsales@mgc-lighting.com

Can't find what you need?
Just call!

Manufacturer	System	Model	Spec	Property	Original or Non Original	MGC Part No.	
Physio Control	Physio Control Defib	Lifepak 15 21330-001176	11.1V 5.7Ah	LiOn	Original	BATT/110344	
Physio Control	Physio Control Defib	Lifepak 9	16V 2.50Ah	Pb	Non original	BATT/110066	
Physio Control	Physio Control Defib	Lifepak 9, 9P- 803704-03	16V 2.50Ah	Pb	Non original	BATT/110095	
Physio Control	Physio Control Defib	Lifepak 12	12V 1.6Ah	NC	Original	BATT/110102	
Physio Control	Physio Control Defib	Lifepak 12	12V 2.50Ah	Pb	Original	BATT/110104	
Physio Control	Physio Control Defib	Lifepak 500, OEM Number 300-5380-030 (not rechargeable)	12V 7.50Ah	Lithium	Original	BATT/110101	
Physio Control	Physio Control Defib	Lifepak 500, OEM Number 300-5380-030 (not rechargeable)	12V 7.50Ah	Lithium	Non original	BATT/110101-A	
Physio Control	Physio Control Defib	Lifepak 1000	12V 4.5Ah	Lithium	Original	BATT/110305-0	
Protocol	Protocol Propaq CS VSM	101-106, 242, 244, 246	8V 2.5Ah	Pb	Non original	BATT/110265	
Welch Allyn	Welch Allyn Vital Signs Monitor	5300/53NOP/53NTO/53OTO 53000/53NTP/5350P	6V 4.5Ah	Pb	Non original	BATT/110921-A	
Welch Allyn	Welch Allyn Vital Signs Monitor	/53NTO/53OTO/53000 53NTP/5350P	6V 4.5Ah	Pb	Original	BATT/110921	

01473 466 300



MGC Lamps Ltd
 1 Sovereign Centre,
 Farthing Road, Ipswich
 Suffolk, IP1 5AP

Manufacturer	System	Model	Spec	Property	Original or Non Original	MGC Part No.	
Welch Allyn	Welch Allyn Vital Signs Monitor	LXI 400732 Type 4500-84	6V 5.5Ah	Pb	Non original	BATT/110922-A	
Welch Allyn	Welch Allyn Vital Signs Monitor	300 Series, Spot 420, 42NOB, 53NTB, 53NTL, 63NTB, 4200-84	6V 4.5Ah	Pb	Non original	BATT/110923-A	
Welch Allyn	Welch Allyn Vital Signs Monitor	Connex VSM6000 (long life version - BATT-99)	11.1V 6.0Ah	LiOn	Original	BATT/110597-0	
Welch Allyn	Welch Allyn Ophthalmoscope Handle	Type 72200, for handle 71000, 71670, 71020	3.5V 0.75Ah	NC	Original	BATT/110910-0	
Welch Allyn	Welch Allyn Ophthalmoscope Handle	Suitable for 72200, 71000, 71670, 71020 handles	3.5V 0.7A	NiMH	Non original	BATT/110910-AN	
Welch Allyn	Welch Allyn handle	Type 72300	3.5V 0.75Ah	NC	Original	BATT/110912	
Welch Allyn	Welch Allyn Spot Vital Signs Monitor	LX1, Type 105632	6.4V 6.0Ah	LiOn	Original	BATT/110922-L	
Zoll	Zoll Defib	NTP2, PD1400/PD1600/PD1700/PD2000/4410/AED PRO Type 8000-0299-10/8000-0299-01	10V 2.50Ah	Pb	Non original	BATT/110087	
Zoll	Zoll Defib	R Series Type Surepower 8019-0535-01	10.8V 5.8Ah	LiOn	Original	BATT/110467-0	

United Kingdom Sales

☎ 01473 466 300
☎ FREEPHONE 0800 137 739
✉ enquiry@mgc-lighting.com






International Sales

☎ +44(0) 1473 462321
✉ internationalsales@mgc-lighting.com

Can't find what you need?
Just call!

Duracell Industrial Batteries



Size	Property	Voltage	Pack Size	MGC Part No.	
AAA (LR03)	Alkaline-Manganese Dioxide	1.5 V	10 per pack	BAT/AAA/ALK/10P:DUR	
AA (LR6)	Alkaline-Manganese Dioxide	1.5 V	10 per pack	BAT/AA/ALK/10P:DUR	
C (LR14)	Alkaline-Manganese Dioxide	1.5 V	10 per pack	BAT/C/ALK/10P:DUR	
D (LR20)	Alkaline-Manganese Dioxide	1.5 V	10 per pack	BAT/D/ALK/10P:DUR	
9V (6LF22)	Alkaline-Manganese Dioxide	9 V	10 per pack	BAT/9V/ALK/10P:DUR	

01473 466 300



MGC Lamps Ltd
 1 Sovereign Centre,
 Farthing Road, Ipswich
 Suffolk, IP1 5AP

Terms and Conditions

We are also able to offer a build or refurbishment service for obsolete or non-standard equipment.

Battery Life

Guidelines for battery shelf life are based on storing in a dry area at an ambient temperature of 20°C.

Chemistry	Shelf Life
Nickel Cadmium	5 years when stored in the open circuit condition
Nickel Metal Hydride	6 months to top up charge - Do not allow to become fully discharged
Sealed Lead Acid	6 months to top up charge - Do not allow to become fully discharged
Alkaline	5 years
Lithium Manganese Dioxide	5 years
Lithium Poly Carbonmonofluoride	10 years
Lithium Thionyl Chloride	10 years
Silver Oxide	2 years
Zinc Air	5 years in unopened packets
Zinc Carbon	2 years
Zinc Chloride	2 years

Warranty

All of our products carry a 12 month warranty from the date of purchase.

Due to the various products we supply, they will be covered as follows:

- * Nicad, Nimh and Li-Ion Batteries shall maintain a minimum of 80% capacity
- * Sealed Lead Acid Batteries shall be free of defects in material and workmanship for 12 months
- * Electronic items carry a 12 month parts and labour warranty
- * Non rechargeable batteries shall perform at 100% rated capacity upon delivery

No returns are accepted without a RMA number, all items returned to MGC under warranty should be accompanied with a short report detailing the nature of the fault. Please ensure you include your name, address and telephone number with the item.

For a returns authorisation, please contact Customer Services on 01473 466300

Other MGC Brochures available:



MGC Fibre Optic Light Guide Leads, Headlights and Testers



MGC Medical Catalogue



MGC Keeler Brochure

United Kingdom Sales

☎ 01473 466 300

☎ FREEPHONE 0800 137 739

✉ enquiry@mgc-lighting.com

International Sales

☎ +44(0) 1473 462321

✉ internationalsales@mgc-lighting.com

Can't find what you need?

Just call!



MGC
MEDICAL

MGC Lamps Ltd
✉ 1 Sovereign Centre, Farthing Road
Ipswich, Suffolk, IP1 5AP

Copyright © 2016/17 MGC Lamps Ltd

40
YEARS OF LIGHTING

United Kingdom Sales
☎ 01473 466 300
☎ FREEPHONE 0800 137 739
✉ enquiry@mgc-lighting.com

International Sales
☎ +44(0) 1473 462321
✉ internationalsales@mgc-lighting.com

